

4

3

2




1

B

B

A


A



NOTE:
• USE "VERTREL XSI" SOLVENT TO EXPAND SILICONE TUBING PRIOR TO INSERTING COILS.



THE INFORMATION CONTAINED IN THIS DRAWING IS PUBLICLY RELEASED BY COSMIC UNDER AN OPEN SOURCE LICENSE. COSMIC IS NOT LIABLE FOR ANY ASPECT OF USE OF THIS DRAWING, INCLUDING EDITS, FABRICATION, AND REPRODUCTION OF THIS DESIGN.



MATERIAL	FINISH
----------	--------

DOC. NAME
SA, Cable Body, Insulated DFT Filars, Blue-Clear

PROJECT NNP	DATE 6/27/2024	SIZE B	DOC. NUMBER NNP-DWG-145-007-001	REV 1.0
SHEET 1 OF 1	SCALE 2:1			

ITEM NO.	DESCRIPTION
1	Filar, DFT, Blue-Clear
2	Cable, Insulating Silicone Tubing
3	Cable Body, Insulated DFT Filars, Blue-Clear

4

3

2

1